# Mg-FixationScrew/Mg-FixationScrew<sup>c</sup> Gebrauchsanweisung/Instructions for Use

Version 003, gedruckt/printed: 10-2020

Meotec GmbH  
Philipsstr. 8, DE-52068 Aachen

Internet / Web: www.meotec.eu  
E-Mail: info@meotec.eu  
Telefon / Phone: +49 241 475 718 0  
Fax: +49 241 475 178 90

## DEUTSCH

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## ENGLISH

It is necessary that all statements and instructions of this user manual are followed and fulfilled.

### Symbols

Manufacturer

Article number

Batch code

Do not use product after

Do not re-use

Do not use if packaging is damaged

Keep away from moisture

Keep away from sunlight

Designation

Sterilized using irradiation

Do not re-sterilize

Consult accompanying instructions for use

MR safe

## 1. Product Overview

Mg-FixationScrew implants, made of a magnesium alloy (Mg-Y-RE-Zr), have metallic properties, gradually degrade within the body, and are replaced by endogenous tissue. The biomechanical properties are very similar to those of human bone.

This instruction for use applies to the product variants of the cortical bone screws (CBS), pins (PIN) and tenting screws (TNS) of the Mg-FixationScrew implants. The Mg-FixationScrew<sup>®</sup> also has a highly innovative ceramic surface finish. The biodegradation of the implant is therefore delayed postoperatively within the first months. The ceramized implants are identified with a superscript "C" in combination with the product name.

## 2. Intended Purpose

The Mg-FixationScrew is a bioabsorbable bone screw that is used to restore the bone continuity by the fixation of bone fragments or tissue substitutes.

### Indications

#### CBS/TNS:

- Intraoral fixation of bone fragments or substitutes to the jawbone,
- intraoral fixation of membranes to the jawbone,
- extraoral fixation of bone fragments or substitutes to the jawbone.
- stable mandibular or maxillary reconstruction under minor load,

- intra- or extraarticular fractures of small bones and bone fragments,
- arthrodosis, osteotomies or pseudarthrosis of small bones and joints,
- small bony ligament and tendon ruptures,
- osteochondral fractures and dissections,
- osteochondrosis dissecans,
- phalangeal and metacarpal bones.

#### PIN:

- Intraoral fixation of bone fragments or substitutes to the jawbone,
- intraoral fixation of membranes to the jawbone,
- extraoral fixation of bone fragments or substitutes to the jawbone.

### Contraindications

The use of Mg-FixationScrew implants is prohibited in specific clinical situations (absolute contraindication) or should only be planned after special consideration (relative contraindication):

#### a. Absolute contraindications

- Insufficient bone mass for anchorage of the implant,
- confirmation or suspected septic infectious surgical site,
- known allergies and/or known foreign body reaction,
- load-bearing stable osteosynthesis.

#### b. Relative contraindications

- Acute sepsis,
- anamnestic evidence of infections or frequent falling incidents,
- insufficient bone quality (e.g., osteoporosis) and/or local bone tumour,
- alcohol, nicotine and/or drug abuse,
- epilepsy,
- poor skin/soft tissue conditions,
- uncooperative patient and/or patient with restricted intellectual capacity.

### 3. Adverse Reactions

In isolated cases, temporary tissue reactions due to the formation of hydrogen gas may be observed after the procedure. In very rare cases, there are hydrogen gas accumulations that are clinically harmless and recede with time.

### 4. Warnings

When using other makes of implants at the same time, it is important to note that steel, titanium and cobalt-chromium alloys in the surgical site must not be in direct contact with Mg-FixationScrew implants for an extended period (physical contact between implants).

To minimize the risk of swallowing Mg-FixationScrew implants, it must first be ensured that the implant fits securely on the screwdriver blade prior to insertion into the patient's oral cavity.

The Mg-FixationScrew implants are intended for single use only. Re-use constitutes gross negligence. Restenilization of the implants will have an unpredictable impact on the product.

### 5. Patient Information

The surgeon must explain to the patient the factors which may impair the success of the operation and the possible adverse effects as described in the sections 'Contraindications' and 'Adverse reactions', insofar as they are applicable to the patient. The patient must also be informed about the measures he/she can take in order to reduce the likelihood of impairments and adverse side effects.

### 6. Directions for Use by the Surgeon

Before the operation, the surgeon must inform himself exactly about the use of the product and the product itself as well as the specifics of the surgical technique. Before use, the suitability of the product for the planned intervention must be checked. Further product-specific information is available from Metec GmbH, Germany.

## 7. Handling

The expiration date must be checked prior to using the Mg-FixationScrew implants. The implants must not be used after the expiration date. Implants that have been removed from the sterile packaging and were not used, must not be resterilized, but discarded. For implant lengths longer than 5 mm, a pre-drilling with a diameter equal to or not more than 0.2 mm smaller than the thread diameter is mandatory. In case of implant lengths smaller than or equal to 5 mm (CBS, PIN) depending on the patient's bone rigidity, a pre-drilling equal to or not more than 0.4 mm smaller than the thread or shaft diameter is recommended. For further details about the cleaning and sterilization procedures for the instrumentation it is to be referred to the respective instructions for Use of the manufacturer.

## 8. Traceability

To ensure full traceability, the self-adhesive identification labels enclosed within the original packaging must be kept and attached to the operation report as well as the patient dossier. In order to assure immediate action in the event of changes or recalls, the person responsible for the hospital's material store must keep a record of the devices in stock, including the article no. (REF) and batch no. (LOT).

## 9. Implant Removal

If in special cases the implants must be removed (e.g., infection), the same instruments used for implantation and additional forceps may be used. In case of a well-advanced degradation, the implants can be drilled out with a conventional drill bit. Special instruments are not required.

## 10. MRI Safety Information

 The Mg-FixationScrew product family is MR safe.

## 11. Packaging & Sterility

The implants are packaged in sterile protective packaging and must be stored in the same packages. The packaging consists of the storage package as well as a sterile-barrier-system. Products marked 'STERILE/R' have been exposed to gamma radiation of 25 kGy (2.5 Mrad). Only intact sterile packaging guarantees the sterility of the Mg-FixationScrew implants. The sterile packaging must not be opened until immediately before use of the implant. The package must be checked for damages prior to opening. The relevant regulations concerning handling of sterile implants have to be complied with when removing or using the implants.

## 12. Storage Conditions

The sterile package and implant must not be exposed to damage by extreme temperature or mechanical influences during storage. The Mg-FixationScrew implants must be stored in their original packaging in clean and dry conditions protected from vermin, extreme temperatures and moisture. The storage conditions for products marked 'STERILE/R' are specified on the packaging label.

## 13. Directive Compliance

This medical device bears the CE-mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE-mark indicates the responsible notified body.

## 14. Disposal

To dispose these medical devices, no special measures are necessary. National laws and regulations must be observed.